


Transition to GMP, chemically-defined, xeno-free cryopreservation media increases post-thaw functionality of clinically-relevant cell banks

Brian J. Hawkins^{1,3}, Alireza Abazari¹, Faith Kalucki², Aby J. Mathew¹, Benjamin Fryer²

¹BioLife Solutions, Bothell, WA, ²Heart Regeneration Program, University of Washington, Seattle, WA, ³Department of Anesthesiology and Pain Medicine, University of Washington, Seattle, WA

BioLife SOLUTIONS
BIOPRESERVATION TOOLS FOR CELLS, TISSUES, AND ORGANS



ISCRM Institute for Stem Cell & Regenerative Medicine
at the University of Washington

ABSTRACT

Cryopreservation is an essential element for the generation of master and working banks for cell therapy and regenerative medicine applications. As these products are scaled-up for commercialization, cryopreservation protocols established early in preclinical development must be amended to comply with Good Manufacturing Practices that allow for large-scale production of cells for clinical use. Optimization of cryopreservation protocols for master and working cell banks prior to scale-up therefore ensures GMP compliance, acceptable viability, and a return-to-function of frozen stocks post-thaw. Using representative human embryonic stem cells (hESC) and Jurkat T-cells, our group demonstrates that utilization of a GMP-manufactured, serum-free, protein-free, animal-origin-free, cryopreservation media (CryoStor[®]) resulted in a significant increase in viable cell recovery and functional cell recovery as determined using dye exclusion assays and functional indicators for up to 48 hours post-thaw. The intracellular-like media also reduced the DMSO concentration required to achieve acceptable post-thaw viable recovery. The process outlined in these studies provides a roadmap that can be rapidly and easily implemented to cryopreserve valuable cell therapy samples, and ensure desired viability and functional recovery of master and working cell banks following long-term frozen storage.

METHODS

Cell Culture: The Jurkat (Clone E6-1) human acute T-cell leukemia (ATCC, Manassas, VA) was cultured in RPMI 1640 (Lonza, Walkersville, MD) supplemented with 10% FBS (Atlas Biologicals, Fort Collins, CO). The ELF1 human Embryonic Stem Cell line (University of Washington, Seattle, WA) was cultured in an albumin-free, chemically defined, E8[™] medium (ThermoFisher Scientific, Waltham, MA).

Cryopreservation: Jurkat and ELF1 cells were cryopreserved in the indicated cryomedia, placed in 1.2 or 2.0 ml cryovials, and cooled at 2-8°C for 10 min. Cryovials were frozen using either an isopropanol freezing device or controlled rate freezer, and then transferred to LN₂ storage for a minimum of 24 hr. Samples were thawed in a 37°C water bath, immediately resuspended in culture media (1:10 dilution), and then either plated on laminin-coated culture vessels (hESCs) or cultured in suspension (Jurkat).

Viability & Recovery: Structural cell viability was evaluated immediately post-thaw and at indicated time periods via membrane integrity using either the NucleoCounter NC-200 or NC-3000 imaging cytometer (ChemoMetec) and a Via-1 cassette.

Apoptosis: Post-cryopreservation caspase 3/7 activity was determined in Jurkat T-cells at the indicated time point using the NucView488[™] fluorescent substrate and the NC-3000 equipped with the Flexicyte package. Results were normalized to non-frozen controls.

Cellular Glutathione: Cellular thiol content (glutathione) was determined using the VitaBright-48 bimane reagent and the NucleoCounter NC-3000.

Statistical Analysis: For all measurements, data represent the mean of 3-14 independent experiments and are normalized to pre-freeze conditions and 37°C controls as indicated. Groups were analyzed for statistical significance via two-way ANOVA with Tukey post-hoc comparisons and p<0.05 set *a priori*. Error bars represent standard error of the mean (SEM).

Abbreviations: Ctrl – culture control, CS5 – CryoStor[®] CS5, PLA – Plasmalyte A, HSA – human serum albumin, GSH – glutathione, DMSO – dimethylsulfoxide, GMP – Good Manufacturing Practices

CONCLUSIONS

- Measures of immediate post-thaw membrane integrity can vary over time and may not accurately reflect long-term cell viability following cryostorage. Additional measures of cellular function should be employed to more accurately assess cell health following cryopreservation.
- The onset of cryopreservation-induced apoptosis can be accurately detected within hours post-thaw by the measurement of caspase 3/7 activity and DNA fragmentation. The presence of apoptotic indicators corresponds to the observed decline in long-term cell viability as determined by membrane integrity.
- Altered cellular redox balance can be detected rapidly post-thaw and suggest oxidative stress following cryopreservation. Oxidative stress and decreased GSH levels may contribute to cryopreservation-induced cell death/damage.
- Measures of post-thaw cell viability and functionality are highly dependent upon the cryomedia employed. Optimized CryoStor[®] cryomedia improves viable recovery versus traditional 'home-brew' cryomedia, reduces cryopreservation-induced apoptosis, and conveys more accuracy to membrane integrity as an indicator of cell viability.
- Optimized CryoStor[®] media exhibits increased viable recovery and functionality post-thaw versus traditional 'home-brew' cryomedia. The effectiveness of CryoStor[®] CS5 was superior to that of a non-GMP traditional 'home-brew' T-cell cryomedia composed of Plasmalyte A/10% w/v Human Serum Albumin/5% v.v DMSO. Using CryoStor[®] media, Jurkat T-cells and ELF1 hESCs can be successfully cryopreserved following a traditional freezing protocol.
- These data stress the need for multiple measures of cell viability post-thaw, and suggest that optimized CryoStor[®] media offers a more robust solution for the accurate evaluation of post-thaw viability from master and working cell banks using multiple measures of viability post-thaw.
- hESCs can be expanded, cryopreserved, thawed and successfully returned to culture in a cell culture media and GMP-manufactured cryomedia that are chemically defined and do not contain serum or albumin. These conditions will be required to manufacture cellular therapies derived from pluripotent cells.
- Intracellular-like, chemically defined, and xeno-free CryoStor[®] media manufactured under GMP can facilitate the rapid optimization of cryopreservation protocols for clinically-relevant master and working cell banks.

RESULTS




Figure 1. CryoStor[®] reduces post-cryopreservation cell loss and speeds cell expansion. Jurkat cells cryopreserved at a density of 5x10⁶/ml in either a traditional T-cell cryomedia composed of Plasmalyte A supplemented with 10% human serum albumin and 5% DMSO (HSA/PLA) or GMP-manufactured CryoStor[®] media preformulated with 5% DMSO (CS5) were thawed in a 37°C water bath after LN₂ storage. Immediately post-thaw and at increasing intervals up to 48 hr, viability and cell number was determined using the Via-1 cassette and visualized by an automated NucleoCounter NC-3000 imaging cytometer. (A) Initial cell recovery/count was similar between CS5 and HSA/PLA but increased significantly in CS5 at 48 hr (**p<0.001). Jurkat cells cryopreserved in CS5 exhibited increased (B) viability (**p<0.001) at 24 hr and (C) viable recovery (****p<0.0001). Jurkat expansion was evident at 48 hr post-thaw following cryopreservation in CS5.




Figure 2. CryoStor[®] reduces cryopreservation-induced apoptosis. To determine the mechanism by which cryopreservation associates with a decline in cell viability at 24 hr, Jurkat cells cryopreserved in the indicated cryomedia were evaluated for caspase 3/7 activity at incremental time points. (A) Increased caspase 3/7 activity was evident 2 hr post-thaw in Jurkat cells cryopreserved in HSA/PLA that normalized following 24 hr of post-thaw culture. Caspase 3/7 activity was similar to non-frozen controls at all time points in Jurkat cells cryopreserved in CS5. Representative image cytometry dot-plot of fluorescence intensity of the caspase 3/7 fluorescence substrate (Nuc488[™]) and the necrotic indicator RedDot in Jurkat cells cryopreserved in (B) CS5 and (C) HSA/PLA. Fluorescence intensity was normalized to non-frozen controls as indicated. (D) Cell cycle analysis using DAPI revealed a significant increase in Jurkat cells exhibiting fragmented DNA when cryopreserved in HSA/PLA versus CS5, and correspond with the downstream effects of caspase 3/7 activity. Dashed lines in (A) and (D) represents caspase 3/7 activity and % cells with fragmented DNA in non-frozen controls, respectively. (****p<0.0001)




Figure 3. CryoStor[®] protects cellular redox balance. Glutathione (GSH) is the largest contributor to the cellular thiol pool and is essential for numerous oxidation/reduction reactions. An oxidative shift in the cellular redox balance results in the dimerization of two glutathione molecules and a general loss of cellular antioxidant capacity. Oxidative stress can be a trigger for apoptosis, secondary necrosis, and cellular dysfunction during environmental stress. To determine whether cryopreservation-induced cell death/damage associates with oxidative stress, Jurkat T-cells were evaluated immediately post-thaw and at increasing time intervals with the thiol-reactive fluorophore VB-48. Fluorescence dot-plots of VB-48 (reduced cellular thiol) versus Propidium Iodide (PI) (plasma membrane permeability) in Jurkat cells cryopreserved in (A) CS5 and (B) HSA/PLA. Fluorescence intensity was normalized to non-frozen controls as indicated. (C) Cryopreservation in HSA/PLA resulted in a large population of cells with low thiol (GSH) content versus non-frozen controls. Cryopreservation in HSA/PLA further required extended post-thaw culture to normalize cellular redox balance. In contrast, CS5 preserved cellular redox balance immediately post-thaw and throughout the 24 hr post-thaw culture.




Figure 4. Optimized cryopreservation of a human Embryonic Stem Cell (hESC) working cell bank for preclinical development. The ELF1 hESC line was cryopreserved in CryoStor[®] cryomedia at increasing DMSO concentrations at a density of ~1x10⁶ cells/ml in the presence of the ROCK1 inhibitor Y-27632. hESCs were frozen using an isopropanol freezing chamber to -80°C and then transferred to liquid phase LN₂ storage. hESCs were thawed in a 37°C water bath and then resuspended at a ratio of 1:9 in CGM. Immediately post-thaw, viability was assessed via membrane permeability using the NucleoCounter NC-100. (A) ELF1 hESCs exhibited elevated cell viability following cryostorage in CryoStor[®] media at all DMSO concentrations, compared to a PBS control, that was maximal using CS10 (10% DMSO). Post-thaw cell culture revealed improved long-term viability using CS10 (data not shown), which was chosen for subsequent cell manufacturing of a GMP working cell bank. hESCs were cryopreserved in CS10 and frozen using a LN₂ controlled rate freezer at a density of 5x10⁶ cells/ml to generate a working cell bank for preclinical development. Following LN₂ storage, hESCs were thawed as described previously and viability assessed using a 21 CFR Part 11-compliant NucleoCounter NC-200. hESCs were plated on laminin-coated culture dishes and expanded for 96 hr before harvest. (B) Microscopic evaluation of hESC morphology at 24 and 96 hr post-thaw reveals enhanced cluster formation and rapid expansion of cells cryopreserved in CS10. (C) Measures of viability, recovery, viability recovery, hESC yield, and fold expansion of hESCs after 96 hr of post-thaw culture. Optimized cryopreservation of hESCs using CS10 resulted in a GMP-compliant working cell bank that can be rapidly expanded for preclinical development.