

# PEACE OF MIND THAT SAMPLES ARE SAFE


SenseAnywhere | Monitor and audit safe storage



**With SenseAnywhere, it's easier than ever to know samples and specimens are safe, whether in the lab or in remote settings.**

SenseAnywhere is a cloud-based tool that allows users to monitor and audit the safe storage of valuable biological specimens remotely — from anywhere in the world — with just an internet connection.

SenseAnywhere combines wireless sensors and the power of the cloud to track ULT temperatures across a few units or an entire fleet.



**Sense**  
**Anywhere**

## Why SenseAnywhere?

- ⦿ Connect and monitor multiple units. Users can connect any number of freezers to view performance data from multiple ULT units from one location on any device.
- ⦿ Maintain compliance. Users will have a record of temperature performance to support Good Manufacturing Practice (GMP) Full 21 CFR Part 11 and EN12830-2018 standards and regulatory compliance.
- ⦿ Secure data. Automatic data uploads (via site Ethernet or cellular GSM) and an independent system for redundancy maintains data safely.
- ⦿ Instant notifications. Users can choose to receive customized text or email alarm notifications and alerts, so they'll know the moment their ULTs dip below temperature thresholds.
- ⦿ Minimize IT roadblocks. SenseAnywhere is easy to set up and rarely requires IT support.
- ⦿ Large temperature display. A large home screen (on any device with a web browser connected to the cloud) makes it easy to see temperature information at a glance. Users can view ambient temperatures, freezer temperatures, temperature history and alarm details.
- ⦿ Monitor other lab equipment. SenseAnywhere doesn't just work with Stirling Ultracold ULTs. Users can monitor temperatures and receive alerts from other lab equipment types and manufacturers.

# STEP 1

## Configure the best SenseAnywhere solution for the unique needs of each user

It's easy to configure the best SenseAnywhere solutions for your customers.

### Identify if a stationary or mobile environment needs to be monitored.

- Stationary monitoring utilizes the Indoor Access Point and mobile requires the Mobile Access Point.

### Layout a floor plan that shows the monitoring locations.

- This helps determine the number of Access Points required.
- Indoor Access Points can communicate up to 600 meters line of sight.
- Doors, walls and liquids may diminish the distance of communication.

## Determine the system components



### Plan Details

Text/SMS	Yes
E-mail/Cloud access	Yes
Temperature data storage period	15 years
Max Sensors viewable in chart exports	100
Max users	Unlimited
Alarm profiles per Sensor	5
Reminder alerts	5
Credits per year, per Sensor	3

# STEP 2


## Set up remote monitoring

Setting up remote temperature monitoring on a Stirling Ultracold ULT freezer is as easy as A-B-C.

# A

### Insert the temperature probe into the freezer\*

Insert the SenseAnywhere PT100 Probe into the freezer. Attach the wireless sensor (AiroSensor X ER) to the side of the freezer near the access port.




# B

### Set up an Indoor or Mobile AccessPoint

**Indoor Access Point:** Plug into AC power and insert the Ethernet cable.

**Mobile Access Point:** Plug into AC or DC power and insert the SIM Card.



# C

### Log into your SenseAnywhere dashboard

Log into the dashboard via any internet-connected device to view real-time metrics and set communication preferences.



\*SenseAnywhere can also be configured for monitoring humidity and optional door opening status with the proper setup. To configure a custom solution and learn more contact a Stirling Ultracold representative.

Or download our [SenseAnywhere Installation Guide](#).



Stirling Ultracold  
6000 Poston Road  
Athens, Ohio 45701, USA  
T +1.740.274.7900 / 855.274.7900  
F +1.740.274.7901  
E [info@stirlingultracold.com](mailto:info@stirlingultracold.com)  
[BioLifeSolutions.com](http://BioLifeSolutions.com)

2023© Stirling Ultracold, a part of BioLife Solutions. All rights reserved.

Global Cooling technology is manufactured under U.S. and International patents. Stirling Ultracold is a trademark of Global Cooling, Inc. Specifications subject to change without notice.

### Storage Solutions

SS-ULT-230630\_R02