

## Best-in-class optimized biopreservation media for cells and tissues

- 100% DMSO
- Sterile Vial
- Standard and Custom Fill Volumes
- USP/Highest Quality Components
- Filled Under cGMP
- Sterility and Endotoxin Testing



BloodStor®, is a series of generic cGMP freeze media products used to cryopreserve stem cells, other cells, and cell components isolated from umbilical cord blood, peripheral blood, bone marrow, and other biologics. All BloodStor® products are formulated using USP grade components and Water for Injection (WFI) quality water.

Test methods and criteria are provided on all lot specific Certificates of Analysis and Release.

#### Glossary of label symbols:

REF	LOT	STERILE
Order Number	Number Batch Code	Aseptic Processing Technique
Exp. Date	Temperature	Manufactured by



3303 Monte Villa Parkway, Suite 310  
Bothell, WA 98021 USA  
1.866.424.6543 Phone | 1.425.402.1433  
BioLifeSolutions.com

### Ordering Information

Product Name	Size	Part #
BloodStor® 100	50mL vial	410301
BloodStor® 100	100mL vial	410302

### Related Products

BloodStor® 55-5  
Please visit website for ordering information.

### Shipping and Storage

BloodStor® 100 ships at ambient temperature. Upon receipt, store at 20° to 30° C, protected from light, until ready to use.

### To Order

**Call:** 1.866.424.6543 | **Fax:** 425.402.1433  
**Sales:** SalesOne@BioLifeSolutions.com  
**Web:** BioLifeSolutions.com  
**Technical Support:** info@BioLifeSolutions.com