

# BloodStor® 55-5

Freeze media

- ✓ 55% (w/v) DMSO
- ✓ 5% (w/v) Dextran 40
- ✓ Standard and Custom Fill Volumes
- ✓ USP/Highest Quality Components
- ✓ cGMP Manufactured
- ✓ Sterility, Endotoxin, and Cell-Based Release Testing



## BEST-IN-CLASS OPTIMIZED BIOPRESERVATION MEDIA FOR CELLS AND TISSUES

**BloodStor®** is a series of generic cGMP freeze media products used to cryopreserve stem cells and other cells isolated from umbilical cord blood, peripheral blood, bone marrow, and other biologics. All BloodStor products are formulated using USP grade components and Water for Injection (WFI) quality water. Standard and custom packaging and fill volumes are available.

### To Order

**Call:** 1.866.424.6543 | **Fax:** 425.402.1433

**Sales:** salesone@BioLifeSolutions.com

**Web:** BioLifeSolutions.com

**Technical Support:** info@BioLifeSolutions.com

### ORDERING INFORMATION

Product Name	Size	Part #
BloodStor® 55-5	5.6mL vial	355371
BloodStor® 55-5	7.2mL vial	355372
BloodStor® 55-5	1L bag	355210
BloodStor® 55-5	8mL syringe	355469

### RELATED PRODUCTS

BloodStor® 100

Please visit website for ordering information.

Materials are manufactured under cGMP

Test methods and criteria are provided on all lot specific Certificates of Analysis and Release.

### SHIPPING AND STORAGE

BloodStor® 55-5 ships at ambient temperature.

Upon receipt, store at 2° to 8°C, protected from light, until ready to use.

6008\_05

3303 Monte Villa Parkway, Suite 310  
Bothell, WA 98021 USA  
1.866.424.6543 Phone | 1.425.402.1433 Fax  
**BioLifeSolutions.com**

 **BioLIFE SOLUTIONS®**  
BIOPRESERVATION TOOLS FOR CELLS, TISSUES, AND ORGANS