

# CELL THAWING MEDIA

10% DEXTRAN 40 (in 0.9% NaCl)

- ✓ Pre-Formulated
- ✓ Serum-Free
- ✓ Protein-Free
- ✓ USP/Highest Quality Components
- ✓ cGMP Manufactured
- ✓ Sterility, Endotoxin, and pH Release Testing



## MEDIA FOR THAWING HUMAN CELLS AND TISSUES FOLLOWING CRYOPRESERVATION

**Cell Thawing Media 10% Dextran 40 (in 0.9% NaCl)** provides Dextran and saline for washing cryopreserved cells and tissues to dilute or remove cryoprotectants.

### ORDERING INFORMATION

Product Name	Size	Part #
Cell Thawing Media 10% Dextran 40 (in 0.9% NaCl)	250mL bag	980203

### RELATED PRODUCTS

Cell Thawing Media 10% Dextran 40 (in 5% Dextrose)

Please visit website for ordering information.

To Order

**Call:** 1.866.424.6543 | **Fax:** 425.402.1433

**Sales:** salesone@BioLifeSolutions.com

**Web:** BioLifeSolutions.com

**Technical Support:** info@BioLifeSolutions.com

6021\_02

# CELL THAWING MEDIA

## Cell Thawing Media Usage Protocol

- Cell Thawing Media ships at ambient temperature.
- Upon receipt, store Cell Thawing Media 10% Dextran 40 (in 0.9% NaCl) at 20° - 25°C until ready to use. Protect from freezing.
- Wipe down outside of container with 70% ethanol before opening.
- Contents are sterile. If seal has been broken, do not use and contact BioLife Solutions.

Materials are manufactured under cGMP

Test methods and criteria are provided on all lot specific Certificates of Analysis and Release.