

# BloodStor<sup>®</sup> 27 NaCl

Freeze media

- ✓ 27% DMSO in Saline
- ✓ Sterile Contents
- ✓ 80mL & 112mL Fill Volume
- ✓ USP/Highest Quality Components
- ✓ cGMP Manufactured
- ✓ Sterility, Endotoxin and DMSO Content Testing



## BEST-IN-CLASS OPTIMIZED BIOPRESERVATION MEDIA FOR CELLS AND TISSUES

**BloodStor<sup>®</sup>** is a series of generic cGMP freeze media products used to cryopreserve stem cells, other cells, and cell components isolated from umbilical cord blood, peripheral blood, bone marrow, and other biologics. All BloodStor<sup>®</sup> products are formulated using USP grade components and Water for Injection (WFI) quality water.

### To Order

**Call:** 1.866.424.6543 | **Fax:** 425.402.1433

**Sales:** salesone@BioLifeSolutions.com

**Web:** BioLifeSolutions.com

**Technical Support:** info@BioLifeSolutions.com

### ORDERING INFORMATION

Product Name	Size	Part #
BloodStor <sup>®</sup> 27 NaCl	80mL bag	327207
BloodStor <sup>®</sup> 27 NaCl	112mL bag	327214

### RELATED PRODUCTS

BloodStor<sup>®</sup> 55-5, BloodStor<sup>®</sup> 100

Please visit website for ordering information.

### Materials are manufactured under cGMP

Test methods and criteria are provided on all lot specific Certificates of Analysis and Release.

### SHIPPING AND STORAGE

BloodStor<sup>®</sup> 27 NaCl ships at ambient temperature.

Upon receipt, store at 2° to 25°C, protected from light, until ready to use.

6023\_02

3303 Monte Villa Parkway, Suite 310  
Bothell, WA 98021 USA  
1.866.424.6543 Phone | 1.425.402.1433 Fax  
**BioLifeSolutions.com**

 **BioLIFE SOLUTIONS<sup>®</sup>**  
BIOPRESERVATION TOOLS FOR CELLS, TISSUES, AND ORGANS