




## Media for thawing human cells and tissues following cryopreservation

- Pre-Formulated
- Serum-Free
- Protein-Free
- USP/Highest Quality Components
- cGMP Manufactured
- Sterility, Endotoxin, and pH Release Testing



Cell Thawing Media 10% Dextran 40 (in 0.9% NaCl) provides Dextran and saline for washing cryopreserved cells and tissues to dilute or remove cryoprotectants.

### Glossary of label symbols:

<b>REF</b>	<b>LOT</b>	<b>STERILE</b>
Order Number	Number Batch Code	Aseptic Processing Technique
	 2°C ↑ 8°C	
Exp. Date	Temperature	Manufactured by

### Ordering Information

Product Name	Size	Part #
Cell Thawing Media 10% Dextran 40 (in 0.09 % NaCl)	250mL bag	980203

### Related Products

Cell Thawing Media 10% Dextran 40(in 5% Dextrose)  
Please visit website for ordering information.



3303 Monte Villa Parkway, Suite 310  
Bothell, WA 98021 USA  
1.866.424.6543 Phone | 1.425.402.1433  
BiolifeSolutions.com

### To Order

**Call:** 1.866.424.6543 | **Fax:** 425.402.1433  
**Sales:** SalesOne@BioLifeSolutions.com  
**Web:** BioLifeSolutions.com  
**Technical Support:** info@BioLifeSolutions.com

# Cell Thawing Media

## Cell Thawing Media Usage Protocol

---

- Cell Thawing Media ships at ambient temperature.
- Upon receipt, store Cell Thawing Media 10% Dextran 40 (in 0.9% NaCl) at 20°-25° C until ready to use. Protect from freezing.
- Wipe down outside of container with 70% ethanol before opening.
- Contents are sterile. If seal has been broken, do not use and contact BioLife Solutions.

---

Materials are manufactured under cGMP

Test methods and criteria are provided on all lot specific Certificates of Analysis and Release.

---